# BIOFLOATTM

# The anti-adhesive surface for spheroid culture



### ADVANTAGES OF SPHEROID CULTURE

- Increased cell contacts
- Pronounced extracellular matrix
- Improved in vitro model



### BIOFLOAT<sup>™</sup> solves your challenges in the area of spheroid cultures

Some challenging spheroid cultures have already been successfully established using the BIOFLOAT<sup>™</sup> cell culture surface (e.g. spheroids from primary hepatocytes). A list of cell lines and cell types successfully tested with BIOFLOAT™ can be found on page 6.

In many areas of biomedical research, in vitro models are indispensable. The most conventional form is the two-dimensional cell culture. Discrepancies often occur when transferring the results to an entire organism. The aim of the three-dimensional cell culture is therefore to close this gap between the in vitro and in vivo situation.

Spheroid cultures offer a simple and cost-effective variant of

3D cell culture. The cells form a three-dimensional cellular aggregate with pronounced cell-cell and cellmatrix contacts.

The new BIOFLOAT cell culture surface makes it possible for you to produce perfect spheroids quickly and reproducibly.

BIOFLOAT is used in a variety of fields such as cancer and stem cell research, in the preclinical phase of drug research and in toxicological studies. Spheroid cultures improve the efficiency and reliability of preclinical cell models in these applications.



Fig. 1: Cells of a fibroblast cell line (373) were seeded in various cell counts on the BIOFLOAT™ cell culture plate. An uncoated plate serves as a control. The results were documented microscopically after three days. It can be clearly seen that spheroids are successfully formed using BIOFLOAT<sup>™</sup>. In addition, the size of the spheroid can be influenced by the number of cells/well. On the uncoated surface however, the fibroblasts can adhere and do not form spheroids.





### WHY BIOFLOAT<sup>™</sup> ?

- Robust coating
- Reliable quality
- Easy handling
- Fast results
- High reproducibility

The polymer coating of the BIOFLOAT<sup>™</sup> surface modifies the plastic surface in a simple way. The inert coating contains molecules that are anchored to the polystyrene surface by strong physical interactions and self-organising capability. This achieves a particularly uniform treatment.

The BIOFLOAT<sup>™</sup> surface features a highly anti-adhesive property. This makes it possible for the cultivated adherent cells to form primarily cell-cell contacts without adhering to the surface of the vessel – it forms an ultra-low attachment surface.

Control

BIOFLOA

Competitor



Spheroids cultivated using the BIOFLOAT<sup>™</sup> surface have a very even round shape. Usually, the formation of exactly one spheroid per well is achieved. Both of these features lead to a high reproducibility of your results. This makes BIOFLOAT<sup>™</sup> particularly suitable for high throughput analyses where it is especially important to examine exactly one symmetrical spheroid per well.

The robustness of the BIOFLOAT<sup>™</sup> coating greatly facilitates daily work. The performance of the BIOFLOAT<sup>™</sup> cell culture surface is not impaired even by multiple washing steps or mechanical action from a pipette tip (see Fig. 2).

Heavily scratched

### BIOFLOAT<sup>™</sup> enables spheroid cultures – rapid, reliable and reproducable

Fig.3: 200 ul of a suspension of 3T3 cells with a concentration of 30,000 cells/ml (corresponding to 6,000 cells/well) was seeded per well. Wells with exactly one spheroid were identified and plotted in percent as a function of incubation time.

The reliable quality of the BIOFLOAT<sup>™</sup> cell culture surface enables the formation of perfect spheroids even for challenging cells. This also includes cells that do not form spheroids on existing products.

circle.



Lightly scratched



### **Rapid spheroid formation**

The BIOFLOAT<sup>™</sup> surface enables rapid spheroid formation. Depending on the cell line or cell type, the formation of the spheroids on the BIOFLOAT<sup>™</sup> surface takes between 2 and 24 hours. Uniform spheroids have been shown to form more rapidly than on most anti-adhesive, cell-repellent surfaces (Fig.3).



Fig.5: 100 ul of a suspension of primary human hepatocytes with a concentration of 25,000 cells/ml (equivalent to 2,500 cells/well) was seeded per well. After spheroid formation, 50ul of medium was exchanged every 48-72 hr.



Day 0 Day 2 Day 4 Day 8

HUMAN HEPATOCYTE SPHEROID CULTURES

Spheroid formation of primary human hepatocytes

### **High reproducibility**

Spheroids formed with the help of the BIOFLOAT<sup>™</sup> surface have a high circularity, which enables a high data consistency (Fig. 4). No deposits, satellite aggregates or irregular shapes are formed, which ensures a high reproducibility.

Fig.4: 200 ul of a suspension of 3T3 cells with a concentration of 30,000 cells/ml (corresponding to 6,000 cells/well) was seeded per well. The relative circularity of the spheroids formed was determined and plotted as a function of time. The higher the value, the rounder the spheroid. A value of 1 would correspond to a perfect



BIOFLOAT<sup>™</sup> leads to successful and reliable spheroid formation even with challenging cells



### The following cells have already been successfully tested for the BIOFLOAT<sup>™</sup>-mediated spheroid culture

The <u>BIOFLOAT<sup>™</sup> plate</u> is individually sterile-packed in an aluminium bag. It is also endotoxin-free and non-cytotoxic.

Name	Description		
3T3	Fibroblasts (M. musculus)		
A431	Squamous cell cancer cell line (H. sapiens)		
B16	Melanoma cell line (M. musculus)		
CaCo-2	Colon carcinoma cell line (H. sapiens, Caucasian)		
Capan-1	Pancreatic adenoma cancer cell line (H. sapiens)		
СНО	Ovarian cell line (C. griseus)		
D492	Epithelial breast cancer cell line (stem cell-like) (H. sapiens)		
D492HER	Tumourigenic breast epithelial stem cell line from D492 cells (H. sapiens)		
DAN-G	Pancreatic cancer cell line (H. sapiens)		
ESCs	Embryonic stem cells (S. scrofa domesticus)		
FAMPAC	Pancreatic adenoma cancer cell line (H. sapiens)		
H1975	Lung adenocarcinoma cell line (H. sapiens)		
H2228	Lung adenocarcinoma cell line (H. sapiens)		
H3122	Lung adenocarcinoma cell line (H. sapiens)		
HCC1433	Breast cancer cell line (H. sapiens)		
HCT-116	Colon cancer cell line (H. sapiens)		
hDPSC	Primary dental pulp stem cells (H. sapiens)		
hDPSC+Panc1	Pancreatic cancer cell line (H. sapiens)		
HEK293	Embryonic kidney cells (H. sapiens)		
HepG2	Hepatoma cell line (H. sapiens)		
HT-29	Colon adenocarcinoma cell line ( <i>H. sapiens</i> , Caucasian)		

Name	Description
huARLT	Immortalised endothelial cells (from HUVEC cells) (H. sapiens)
HuOB	Immortalised osteoblasts (H. sapiens)
HUVEC	Venous endothelial cells (H. sapiens)
iPSC-Gata6	iPSC-derived hepatocytes
MCF10A	Breast cancer cell line (H. sapiens)
MCF-7	Breast cancer cell line (H. sapiens)
MDA-MB231	Breast cancer cell line (H. sapiens)
Mia-Paca	Pancreatic cell line (H. sapiens)
Panc1	Pancreatic cell line (H. sapiens)
Panc39	Pancreatic cell line (H. sapiens)
PRH with RHSteC	Hepatic stellate cells/Ito cells (R. norwegicus)
PRH+ HHSteC	Hepatic stellate cells/Ito cells (H. sapiens)
RPMI	B-lymphocyte cell line from myeloma patients (H. sapiens)
SFFV2	Immortalised astrocytes (H. sapiens)
-	Differentiated fatty cell organoids from pluripotent stem cells
-	Endometrial organoids from detached primary cells (non-human primates)
-	Fibroblast progenitor cells (M. cerebralis)
-	iPSC-derived cardiomyocytes (H. sapiens)
-	Liver organoids (differentiated) (M. musculus)
-	Neuronal stem cells (HN9 differentiated)
-	Primary hepatocytes (H. sapiens, M. musculus, M. fascicularis, C. lupus familiaris)



### Order information

SKU	Name	Number of wells	Bottom shape	Packaging
F202003	BIOFLOAT™ 96-well cell culture plates, surface : ultra-low attachment	<u>96</u>	U	4PCs (Sterile)
F202301	BIOFLOAT™ 96-well cell culture (F bottom) plates, surface : ultra-low attachment	<u>96</u>	U	4PCs (Sterile)
F224384	BIOFLOAT™ 384-well cell culture plates, surface : ultra-low attachment	<u>384</u>	F	1PCs (Sterile)
F202005	BIOFLOAT™ Flex coating solution, : ultra-low attachment	-	-	1PCs (Sterile)
F202302	BIOFLOAT™ 6-well cell culture plates, surface : ultra-low attachment	6	F	1PCs (Sterile)
F202303	BIOFLOAT™ 24-well cell culture plates, surface : ultra-low attachment	24	F	1PCs (Sterile)

## faCellitate Gmbh

Franz-Volhard-Straße 5 68167 Mannheim Germany

info@facellitate.com





### If you have any questions: We'd be happy to help!

visit our website www.facellitate.com

