

## BIOFLOAT™ FLEX coating solution


### Coating protocol

#### Protocol to prepare a BIOFLOAT FLEX coated 96-well plate:

For optimal performance we recommend the use of a non-treated polystyrene surface

- 1) Add 100µl faCellitate BIOFLOAT™ FLEX coating solution per well in a 96-well plate
- 2) Wait for 3 minutes
- 3) Aspirate liquid completely
- 4) Wait for 30 minutes for the plate to air-dry at room temperature, within the laminar flow hood
- 5) Coated plates can be used immediately for cell culture or stored at room temperature for later use
- 6) Seed cells using a cell seeding density according to your assay

▲ Please do not use tissue culture treated labware!



For additional information, please visit [www.faCellitate.com](http://www.faCellitate.com) or contact us at

faCellitate  
a venture of Chemovator GmbH,  
Industriestraße 35  
68169 Mannheim, Germany  
T +49 621 60 29872  
info@facellitate.com