



**faCellitate GmbH**  
Industriestraße 35  
68169 Mannheim  
Germany

Combine Function with Your Labware

Product passed visual inspections. Spheroid formation assay was performed using mouse fibroblast 3T3 cells. The seeding density was 3000 cells/well. Single spheroids formed in >95% of wells of the plates tested.

**Sterility:** Tested for sterility according to Ph. Eur. 2.6.1 & USP.

**Endotoxicity:** Tested and met the established criteria. The acceptance level for product is  $\leq 0.04$  EU/mL or  $\leq 4$  EU/device.

**Visual Attributes:** Pass

**Packaging:** Pass

The product is for single use only. The quality of the product cannot be guaranteed after the exp. date.

The product should be stored unopened under dry conditions at room temperature.

faCellitate GmbH ☎ +49 621 877 5353 0 ✉ info@facellitate.com 🌐 <http://www.facellitate.com>  
VAT: DE346773022

Registration Court: Amtsgericht Mannheim, HRB 741442

Managing Director: Dr. Simon Widmaier, Agnes Sackarendt

Bank: Commerzbank, IBAN: DE04545400330205330400, SWIFT: COBADEFFXXX